



MetAmino®

Sample Preparation Kit



 **CHROMSERVIS**

MetAmino® kits offer an easy sample preparation method for your LC-MS or GC-MS analysis. MetAmino® kits include derivatization reagents and all clean-up accessories to prepare your sample for injection. They eliminate time consuming sample preparation procedures. The new clean-up procedure uses a special material as the end-step including filtration (no need to use a syringe filter prior injection). The other advantage is that the derivatization procedure enables you to **extend the analyte list**. Contact us for further details.

Easy sample preparation
using MSPE

75 amino acids,
polyamines, biogen
amines and coenzymes
in 7 to 12 minutes



FAST



EASY



CLEAN

Unique clean-up step

LC-MS
and GC-MS
kit



COMPACT




CLEAR


NIST library
for GC-MS
available

LC-MS separation


The chromatogram shown below displays TIC and some of the analyte's MRMs.




Chromatogram of 24 amino acids




2-ABA (259.7–260.7)




PRO (271.7–272.7)




MET (305.7–306.7)




ORN (388.7–389.7)




TRP (360.7–361.7)




HIS (411.7–412.7)



LYS (420.7–403.7)



PHE (321.7–322.7)



ASP (345.7–346.7)

METAMINO®

METAMINO®


LC-MS kit

| No. | Name | Synonyms | Quantification | [min] | [M+H] ⁺ | Precursor | Quantifier | CE | Qualifier | CE |
|-----|----------------------------|-------------------|----------------|-------|--------------------|-----------|------------|----|--------------|--------|
| 1 | Cotinine | Syringe standard | | 1.15 | 177.1024 | 177.1 | 98 | 24 | 80 | 28 |
| 2 | Putrescine | | | 1.40 | 189.1599 | 189.2 | 72.1 | 16 | 115.9 | 12 |
| 3 | Cadaverine | | | 1.48 | 203.1755 | 203.2 | 86 | 20 | 69.2 | 28 |
| 4 | Homoserine | HSER | | 1.67 | 202.1074 | 202.1 | 146 | 4 | 102 | 8 |
| 5 | Pyroglutamic acid | PGLU | | 1.86 | 186.1125 | 186.1 | 84 | 20 | 130 | 8 |
| 6 | Arginine | ARG | ✓ | 3.35 | 331.2341 | 331.2 | 70.1 | 48 | 115.9 | 32 |
| 7 | Homoarginine | Internal standard | | 3.63 | 345.2497 | 345.2 | 84 | 52 | 125.9 | 24 |
| 8 | Glutamine | GLN | ✓ | 4.12 | 303.1916 | 303.2 | 186 | 12 | 84 | 44 |
| 9 | Anserine | ANS | | 4.20 | 397.2447 | 397.2 | 109.1 | 52 | 226.2 | 28 |
| 10 | Citrulline | CIT | ✓ | 4.20 | 332.2181 | 332.2 | 70.2 | 36 | 215 | 16 |
| 11 | Methionine sulfoxide | MET-SO | | 4.29 | 322.1683 | 322.2 | 248 | 12 | 100 | 32 |
| 12 | Methionine sulfone | MET-SO2 | | 4.30 | 338.1632 | 338.2 | 238 | 12 | 182 | 16 |
| 13 | 5-Hydroxylysine | HLY | ✓ | 4.40 | 345.2021 | 345.2 | 82 | 48 | 128.1 | 28 |
| 14 | 1-Methylhistidine | 1MHIS | ✓ | 4.52 | 326.2075 | 326.2 | 124 | 36 | 224 | 20 |
| 15 | 3-Methylhistidine | 3MHIS | ✓ | 4.55 | 326.2075 | 326.2 | 95.7 | 52 | 270.1 | 24 |
| 16 | Prolylhydroxyproline | PHP | ✓ | 4.71 | 385.2333 | 385.2 | 170.1 | 16 | 114, 70.2 | 36, 56 |
| 17 | Serine | SER | ✓ | 4.73 | 262.1615 | 262.2 | 106 | 16 | 60.1 | 40 |
| 18 | Asparagine | ASN | ✓ | 4.78 | 271.1649 | 271.2 | 240.1 | 16 | 254.1 | 16 |
| 19 | 4-Hydroxyproline | HYP | ✓ | 5.23 | 288.1807 | 288.2 | 86.1 | 32 | 188 | 12 |
| 20 | Glycine | GLY | ✓ | 5.23 | 232.1544 | 232.2 | 76 | 12 | 132 | 4 |
| 21 | N-Acetylaspartic acid | | | 5.36 | 288.1807 | 288.2 | 88 | 28 | 144.1, 158 | 20, 8 |
| 22 | Glycylproline | GPR | ✓ | 5.39 | 329.2072 | 329.2 | 70.1 | 52 | 172.1, 115.9 | 16, 36 |
| 23 | Threonine | THR | ✓ | 5.45 | 276.1806 | 276.2 | 74.2 | 28 | 176.1 | 12 |
| 24 | 5-Aminolevulinic acid | 5-ALA | | 5.63 | 288.1807 | 288.2 | 158 | 8 | 86.1 | 32 |
| 25 | Ethanolamine | EAM | | 6.09 | 262.1649 | 262.2 | 88 | 16 | 144.1 | 4 |
| 26 | Beta-Alanine | 3-ALA | | 6.22 | 246.1701 | 246.2 | 116 | 12 | 90 | 12 |
| 27 | Alanine | ALA | ✓ | 6.25 | 246.1701 | 246.2 | 90 | 16 | 146 | 8 |
| 28 | Spermine | | | 6.39 | 503.3804 | 503.4 | 102 | 44 | 229.2 | 32 |
| 29 | Histamine | HTA | | 6.52 | 312.1919 | 312.2 | 212 | 12 | 95, 112 | 40, 20 |
| 30 | Indoleacetic acid | IAA | | 6.65 | 232.1333 | 232.1 | 130 | 24 | 176 | 8 |
| 31 | gamma-Aminobutyric acid | GABA | ✓ | 6.75 | 260.1857 | 260.2 | 87.1 | 24 | 86 | 24 |
| 32 | Sarcosine | SAR | ✓ | 6.75 | 246.1701 | 246.2 | 90 | 12 | 146.1 | 8 |
| 33 | 4-Aminobenzoic acid | PABA | | 7.00 | 294.1701 | 294.2 | 134.1 | 20 | 91 | 56 |
| 34 | β-aminoisobutyric acid | BAIBA | ✓ | 7.05 | 260.1857 | 260.2 | 130 | 12 | 112 | 20 |
| 35 | 2-aminobutyric acid | 2-ABA | ✓ | 7.35 | 260.1857 | 260.2 | 160.2 | 8 | 104 | 12 |
| 36 | Proline | PRO | ✓ | 7.58 | 272.1857 | 272.2 | 70.1 | 36 | 172 | 12 |
| 37 | Methionine | MET | ✓ | 7.65 | 306.1735 | 306.2 | 204 | 8 | 104 | 20 |
| 38 | Methionine-d3 | Internal standard | | 7.65 | 309.1922 | 309.2 | 207 | 8 | 107 | 16 |
| 39 | Thiaproline | TPR | ✓ | 8.03 | 290.1422 | 290.1 | 88 | 28 | 134, 190.1 | 16, 8 |
| 40 | Asparatame | | | 8.15 | 451.2439 | 451.2 | 120 | 48 | 88 | 44 |
| 41 | Serotonine | | | 8.22 | 377.2072 | 377.2 | 160 | 32 | 303.1 | 8 |
| 42 | 2,4-diaminobutyric acid | DABA | | 8.25 | 375.2490 | 375.2 | 201.3 | 16 | 245 | 12 |
| 43 | Valine | VAL | ✓ | 8.30 | 274.2014 | 274.2 | 72 | 32 | 116 | 12 |
| 44 | Norvaline | | | 8.35 | 274.2014 | 274.2 | 72.1 | 24 | 174 | 8 |
| 45 | Alanyl-lysine | ALA-LYS | | 8.36 | 474.3174 | 474.3 | 84.1 | 60 | 400.2 | 12 |
| 46 | Carnosine | CAR | | 8.36 | 483.2814 | 483.3 | 110.1 | 40 | 212.1 | 24 |
| 47 | Ornithine | ORN | ✓ | 8.45 | 389.2647 | 389.3 | 70.1 | 60 | 315.2 | 8 |
| 48 | Tryptophan | TRP | ✓ | 8.45 | 361.2123 | 361.2 | 259 | 12 | 159 | 28 |
| 49 | Ethionine | ETH | | 8.60 | 320.1890 | 320.2 | 218.1 | 8 | 75.1 | 44 |
| 50 | Histidine | HIS | ✓ | 8.90 | 412.2443 | 412.2 | 110 | 44 | 312.1 | 16 |
| 51 | Lysine | LYS | ✓ | 8.92 | 403.2804 | 403.3 | 84.1 | 56 | 329.1 | 8 |
| 52 | Phenylalanine | PHE | ✓ | 9.04 | 322.2014 | 322.2 | 120 | 36 | 164 | 16 |
| 53 | Leucine | LEU | ✓ | 9.06 | 288.2170 | 288.2 | 86.1 | 20 | 188.2, 130.1 | 8, 12 |
| 54 | Aspartic acid | ASP | ✓ | 9.06 | 346.2225 | 346.2 | 88.1 | 24 | 159.9 | 12 |
| 55 | Spermidine | | | 9.11 | 446.3226 | 446.3 | 198.1 | 28 | 298.1 | 12 |
| 56 | Glutamic acid | GLU | ✓ | 9.17 | 360.2382 | 360.2 | 163 | 16 | 105 | 36 |
| 57 | Allo-isoleucine | aiLE | | 9.22 | 288.2169 | 288.2 | 130.1 | 12 | 86.1 | 24 |
| 58 | Isoleucine | iLE | ✓ | 9.24 | 288.2170 | 288.2 | 130.1 | 16 | 86.1 | 24 |
| 59 | Norleucine | NLEU | | 9.40 | 288.2170 | 288.2 | 86.2 | 24 | 130.1 | 12 |
| 60 | Pipecolic acid | PIP | | 9.45 | 286.2013 | 286.2 | 128.1 | 12 | 84.1 | 40 |
| 61 | Homophenylalanine | Internal standard | | 9.67 | 336.2170 | 336.2 | 91 | 56 | 117 | 32 |
| 62 | 2-Amino adipic acid | AAA | ✓ | 9.67 | 374.2538 | 374.3 | 98 | 32 | 172 | 16 |
| 63 | Adrenaline | ADN | | 10.15 | 484.2541 | 484.3 | 166 | 36 | 466.2 | 4 |
| 64 | Cysteine | CYS | | 10.16 | 378.1946 | 378.2 | 120.1 | 32 | 278 | 12 |
| 65 | 2-Aminoheptanedioic acid | APA | ✓ | 10.20 | 388.2695 | 388.3 | 112 | 36 | 186.1 | 20 |
| 66 | Glutathione | GSH | | 10.27 | 620.3213 | 620.3 | 186.1 | 28 | 130.2 | 44 |
| 67 | Dopamine | DAM | | 10.27 | 454.2435 | 454.2 | 196.1 | 24 | 152 | 40 |
| 68 | Glutamyl-lysine | GLU-LYS | | 10.44 | 588.3855 | 588.4 | 84.1 | 56 | 128.1, 157.9 | 52, 32 |
| 69 | Homocysteine | HCYS | | 10.48 | 392.2103 | 392.2 | 292.2 | 12 | 118.1 | 16 |
| 70 | Diaminopimelic acid | | | 10.54 | 503.3328 | 503.3 | 82.1 | 60 | 127.9 | 48 |
| 71 | Tyrosine | TYR | ✓ | 10.61 | 438.2487 | 438.2 | 136 | 44 | 179.9 | 28 |
| 72 | Cystathionine | CTH | ✓ | 10.68 | 535.3049 | 535.3 | 201.9 | 24 | 88.1 | 56 |
| 73 | Cystine | C-C | ✓ | 10.76 | 553.2613 | 553.3 | 276.1 | 16 | 219.9 | 24 |
| 74 | Kynurenic acid | | | 10.78 | 346.1649 | 346.2 | 246 | 8 | 144, 190 | 52, 24 |
| 75 | Selenocystine | Se-C-C | | 11.12 | 649.1501 | 649.2 | 323.9 | 20 | 575 | 12 |
| 76 | Kynurenine | | | 11.17 | 465.2596 | 465.3 | 146 | 44 | 274.2 | 12 |
| 77 | Homocystine | HC-CH | | 11.28 | 581.2926 | 581.3 | 290 | 16 | 190 | 20 |
| 78 | 3,4-dihydroxyphenylalanine | DOPA | | 11.55 | 554.2961 | 554.3 | 152.1 | 44 | 454.1 | 12 |
| 79 | Theanine | THE | | | 331.2228 | 331.2 | On request | | | |


Note: CE = Collision energy

GC-MS separation


A typical chromatogram of the amino acid standard solution included in this kit.




Chromatogram of standard solution (above)




LEU (EIC 312)




ILE (EIC 283)




PRO (EIC 296)




ASP (EIC 254)




GLU (EIC 282)




MET (EIC 538)



HYP (EIC 312)



PHE (EIC 91)



TRP (EIC 130)

Note: EIC = Extracted Ion Chromatogram.

METAMINO[®]

METAMINO[®]


GC-MS kit

| Peak No. | Name | Synonyms | Quantification | [min] | Mmi* | Diagnostic ions | | |
|----------|-----------------------------------|-------------------|----------------|-------|-----------|-----------------|-----|-----|
| | | | | | | m/z | m/z | m/z |
| 1 | Sarcosine | SAR | | 4.62 | 497.0308 | 270 | 226 | 69 |
| 2 | Alanine | ALA | ✓ | 4.73 | 497.0308 | 270 | 70 | 69 |
| 3 | N-Acetylglycine | | | 4.92 | 299.0392 | 256 | 113 | 212 |
| 4 | Glycine | GLY | ✓ | 4.93 | 483.0151 | 256 | 212 | 56 |
| 5 | 3-Methylglutaric acid | | | 4.97 | 510.0512 | 311 | 282 | 227 |
| 6 | 2-Aminobutyric acid | ABA | | 5.17 | 511.0465 | 284 | 84 | 113 |
| 7 | 3-Hydroxymethylglutaric acid | | | 5.34 | 752.0325 | 285 | 311 | 85 |
| 8 | 3-Alanine | | | 5.36 | 497.0308 | 269 | 256 | 98 |
| 9 | Valine | VAL | ✓ | 5.41 | 525.0621 | 298 | 283 | 98 |
| 10 | N-Acetylcysteine | | | 5.46 | 571.0134 | 282 | 309 | 509 |
| 11 | 3-Aminoisobutyric acid | BAIBA | | 5.47 | 511.0465 | 256 | 112 | 113 |
| 12 | 3-amino-n-butyric acid | | | 5.52 | 511.0465 | 270 | 227 | 283 |
| 13 | Norvalin | NVAL | | 5.69 | 525.0621 | 298 | 256 | 98 |
| 14 | Leucine | LEU | ✓ | 5.92 | 539.0778 | 312 | 270 | 256 |
| 15 | Ethanolamine | EAM | | 5.98 | 513.0257 | 256 | 269 | 270 |
| 16 | Isoleucine | iLE | ✓ | 5.99 | 539.0778 | 283 | 312 | 256 |
| 17 | Allo-isoleucine | aiLE | | 6.02 | 539.0778 | 312 | 283 | 256 |
| 18 | Homoserine | HSER | | 6.04 | 753.0278 | 100 | 283 | 128 |
| 19 | Proline | PRO | ✓ | 6.05 | 523.0465 | 296 | 297 | 69 |
| 20 | 4-Aminobutyric acid | GABA | | 6.20 | 511.0465 | 112 | 256 | 69 |
| 21 | Threonine | THR | ✓ | 6.24 | 527.0415 | 100 | 283 | 483 |
| 22 | Norleucine | NLEU | | 6.28 | 539.0778 | 312 | 256 | 112 |
| 23 | Pipecolic acid | PIP | | 6.41 | 537.0621 | 310 | 518 | 407 |
| 24 | N-Acetylaspatic acid | | | 6.46 | 539.0414 | 270 | 312 | 228 |
| 25 | Asparagine | ASN | | 6.70 | 522.0261 | 295 | 496 | 113 |
| 26 | Aspartic acid | ASP | ✓ | 6.73 | 723.0173 | 254 | 496 | 296 |
| 27 | 3-Methylcysteine | | | 6.76 | 543.0185 | 61 | 300 | 316 |
| 28 | Thioprolin | TPR | | 6.78 | 541.0029 | 314 | 287 | 86 |
| 29 | 2-Hydroxyglutaric acid | | | 6.78 | 738.0169 | 283 | 239 | 511 |
| 30 | Serine | SER | ✓ | 7.07 | 739.0122 | 268 | 295 | 51 |
| 31 | Acetylserine | | | 7.09 | 555.0363 | 268 | 312 | 113 |
| 32 | Pyroglutamic acid | | | 7.43 | 537.0000 | 282 | 310 | 510 |
| 33 | N-Acetylglutamic acid | | | 7.44 | 567.0727 | 282 | 310 | 510 |
| 34 | Glutamic acid | GLU | ✓ | 7.45 | 737.0329 | 282 | 310 | 510 |
| 35 | 4-Hydroxyproline | HYP | ✓ | 7.47 | 539.0414 | 312 | 294 | 68 |
| 36 | Methionine | MET | ✓ | 7.49 | 557.0342 | 538 | 294 | 494 |
| 37 | Cysteine | CYS | ✓ | 7.84 | 754.9895 | 328 | 285 | 113 |
| 38 | Selenomethionine | | | 7.88 | 604.9786 | 282 | 510 | 405 |
| 39 | Ethionine | ETH | | 7.90 | 571.0498 | 282 | 311 | 571 |
| 40 | 2-Amino adipic acid | AAA | | 8.09 | 751.0486 | 124 | 282 | 324 |
| 41 | Phenylalanine | PHE | ✓ | 8.09 | 573.0621 | 91 | 330 | 92 |
| 42 | 3-Hydroxyproline | HYP | | 8.12 | 539.0414 | 312 | 129 | 256 |
| 43 | 2,4-Diaminobutyric acid | DABA | | 8.52 | 752.0438 | 282 | 256 | 325 |
| 44 | 5-carboxymethyl-cysteine | | | 8.55 | 769.0050 | 213 | 259 | 314 |
| 45 | Homocysteine | HCYS | | 8.59 | 769.0045 | 282 | 342 | 82 |
| 46 | 2-Aminopimelic acid | APA | | 8.66 | 765.0642 | 338 | 138 | 95 |
| 47 | Homophenylalanine | Internal standard | | 8.83 | 587.0778 | 283 | 117 | 483 |
| 48 | Histamine | HTA | | 8.78 | 563.0526 | 308 | 320 | 113 |
| 49 | Glutamine | GLN | | 8.86 | 554.0523 | 84 | 282 | 327 |
| 50 | 4-Aminobenzoic acid | PABA | | 9.05 | 545.0308 | 146 | 345 | 346 |
| 51 | 1-Methylhistidine | | | 9.06 | 577.0682 | 95 | 150 | 350 |
| 52 | Chloro-phenylalanine | | | 9.16 | 607.0231 | 125 | 364 | 180 |
| 53 | Methionine sulfone | MET-SO2 | | 9.18 | 589.0240 | 282 | 82 | 189 |
| 54 | Ornithine | ORN | | 9.25 | 766.0595 | 296 | 256 | 69 |
| 55 | Acetyltyrosine | | | 9.49 | 871.0489 | 333 | 289 | 188 |
| 56 | Histidine | HIS | ✓ | 9.55 | 789.0391 | 307 | 362 | 113 |
| 57 | Glycylproline | GPR | | 9.69 | 580.0679 | 70 | 153 | 296 |
| 58 | Lysine | LYS | ✓ | 9.77 | 780.0751 | 310 | 256 | 153 |
| 59 | Tyramine | | | 9.85 | 589.0570 | 346 | 333 | 289 |
| 60 | 2,6-Diaminopimelic acid (isomers) | | | 10.28 | 1006.0616 | 308 | 536 | |
| 61 | Tyrosine | TYR | ✓ | 10.41 | 815.0435 | 333 | 289 | 113 |
| 62 | 5-Hydroxylysine (isomers) | HLY | ✓ | 10.76 | 1022.0565 | 269 | 256 | 69 |
| 63 | Cystathionine | CTH | | 11.12 | 1038.0337 | 328 | 282 | 69 |
| 64 | Dopamine | DAM | | 11.12 | 831.0000 | 256 | 531 | 113 |
| 65 | Tryptophan | TRP | ✓ | 11.19 | 612.0730 | 130 | 131 | 385 |
| 66 | 3,4-Dihydroxyphenylalanine | DOPA | | 11.32 | 1057.0249 | 149 | 531 | 575 |
| 67 | Prolylhydroxyproline | PHP | | 11.52 | 862.0806 | 296 | 297 | 294 |
| 68 | 3-Nitrotyrosine | | | 11.68 | 860.0286 | 334 | 113 | 378 |
| 69 | Selenocystine | Se-C-C | | 11.93 | 1151.879 | 496 | 295 | 268 |

* Mmi = monoisotopic mass of the derivative of the relevant metabolite.

METAMINO[®]

Sample preparation protocol



METAMINO®

METAMINO®

Ordering information

| Sample preparation kit (CF) | Part number | Samples |
|---|----------------|---------|
| MetAmino® sample preparation LC-MS start-up kit, CF | MAK-5857-AA01 | 100 |
| MetAmino® sample preparation GC-MS start-up kit, CF | MAK-5857-BA01* | 100 |
| MetAmino® sample preparation LC-MS basic kit, CF | MAK-5857-CA04 | 400 |
| MetAmino® sample preparation GC-MS basic kit, CF | MAK-5857-DA04* | 400 |
| MetAmino® sample preparation LC-MS reagents kit | MAK-5857-L002 | 100 |

* Coming soon.

Metabolite sample preparation kit (CF)

Derivatization reagents
 Derivatization centrifugal filters
 (set of 100/400 pcs)
 2ml screw top vials (set of 100/400 pcs)
 Micro tube rack
 LC-MS or GC-MS column
 Sample preparation manual

| Sample preparation kit (SPE) | Part number | Samples |
|---|----------------|---------|
| MetAmino® sample preparation LC-MS kit, SPE | MAK-5857-CS01* | 100 |

* Coming soon.

Metabolite sample preparation kit (SPE)

- Derivatization reagents
- Derivatization SPE MetAmino® columns (set of 100 pcs)
- 2ml screw top vials (set of 100 pcs)
- Micro tube rack
- LC-MS column
- Sample preparation manual

Note: The SPE metabolite kit has been developed for users who prefer the SPE technique and do not have the option to use the much easier MetAmino® centrifugal filters.

Recommended accessories



Ohaus centrifuge FC5306



Ohaus Mini Vortex Mixer



MAK-5857-AA01



MAK-5857-CA04



MAK-5857-L002

Benefits

- No need to heat up the sample
- Reagent storage at 4 °C (no freezer required)
- SPE pipetting step replaced by MSPE
- Further extension according to requests from labs
- Ready-to-inject sample after preparation protocol
- Short sample preparation time



Centrifugal tube (OuterTube)
Microspin Filter (Inner Tube)
MetAmino® sorbent bed (inc. 0.22 µm membrane)

Application areas

- Biotechnology
- Clinical/toxicology
- Environmental research
- Food and beverages
- QC of feed production

MSPE – unique technology

- Sample filtration included
- Simplified technique saving your time
- SPE manifold replaced by use of centrifugal filters
- Evaporation step removed

Matrices



Distributor:



Please contact your local distributor for more information.
info@chromservis.eu
www.chromservis.eu

MetAmino® is a trademark registered by Chromservis s.r.o. Other trademarks and trade names may be used in this document to refer to either the entities claiming the marks and/or names or their products and are the property of their respective owners. We disclaim proprietary interest in the marks and names of others.